

INTERFACE

Lode Valiant – GE Medical MAC 5500

The MAC-5500 can only control the GE medical treadmill

For questions and comments, you can contact the Lode service department by e-mail at service@lode.nl or by telephone at +31 50 5712811.

Lode B.V.
Zernikepark 16
9747 AN Groningen
The Netherlands
Tel.: 31(0)50 5712811
Fax: 31(0)50 5716746

ISO 9001:2000 and ISO 13485:2003 certified
e-mail: ask@lode.nl
<http://www.lode.nl>

